

Customer Name \_\_\_\_\_ Customer No. \_\_\_\_\_ Purchase Order No. \_\_\_\_\_ Fax \_\_\_\_\_

Billing Address \_\_\_\_\_ Patient Name \_\_\_\_\_

Shipping Address \_\_\_\_\_ Order Date \_\_\_\_\_ Measured By \_\_\_\_\_

Telephone \_\_\_\_\_ Email \_\_\_\_\_ Shipping Method \_\_\_\_\_

Credit Card Info \_\_\_\_\_

**upper leg kit:** Check the box next to the corresponding length and width.

\*cG cannot be more than 30cm larger than cE2.

	Length			Width		Quantity
	Short	Standard	Long	Regular	Wide	
Left	ℓ D1-G < 37cm	ℓ D1-G 37 - 41cm	ℓ D1-G > 41	cG < 80cm	cG < 125cm	
Right				cE2 < 65cm	cE2 < 110cm	

Each upper leg kit includes:

- 1 reduction upper leg component
- 1 reduction shelf strap
- 6 hook and loop stays
- 1 Built-In-Tension system guide card
- 1 paper measuring tape
- 1 direction for use

**knee kit:** (the reduction knee kit is only available in one length and width)

\*cE2 cannot be more than 30cm larger than cC.

	Length	Width	Quantity
Left	30cm	cE2 < 110cm cE < 102cm cC < 80cm	
Right			

Each knee kit includes:

- 1 reduction knee component
- 2 hook and loop knee spines
- 1 paper measuring tape
- 1 direction for use
- 1 reduction shelf strap

**lower leg kit:** Check the box next to the corresponding length and width.

\*cC cannot be more than 30cm larger than cB.

	Length			Width		Quantity
	Short	Standard	Long	Regular	Wide	
Left	ℓ mall-D1 < 32cm	ℓ mall-D1 32 - 36cm	ℓ mall-D1 > 36cm	cC < 65cm	cC < 80cm	
Right				cB < 55cm	cB < 65cm	


Each lower leg kit includes:

- 1 reduction lower leg component
- 1 pair reduction kit undersleeve leg
- 1 reduction shelf strap
- 6 hook and loop stays
- 1 Built-In-Tension system guide card
- 1 paper measuring tape
- 1 direction for use

**foot options:** Indicate quantity in appropriate check boxes.

\*foot options are not included with kits and must be ordered separately.

	customizable interlocking afw	compression anklet		compressive undersocks		pac band™	single band afw
	One Size	Standard	Large	Standard	Large	X-Large	One Size
cA1	< 44cm	19-36 cm	22-41 cm	19-32 cm	24-41 cm	< 44cm	
cB	< 56cm	19-27 cm	25-42 cm	19-27 cm	25-42 cm	< 56cm	
Left							
Right							



**reduction kit extension**

Short (30 cm)	Standard (35 cm)	Long (40 cm)

**additional ancillary items:**

x-wide undersleeve leg Quantity (Pairs)	reduction kit shelf strap Quantity (Each)

Leg Measurement Guide	
G	Groin
E2	15cm above E
E1	5cm above E
E	Center of patella
D1	Knee crease
D	Slightly below knee
C	Widest part of calf
B	Narrowest part of ankle above malleolus
A1	Middle of foot
A	Ball of foot
mall	Center of malleolus
	ℓ = length in a straight line
	c = circumference

Customer Name \_\_\_\_\_ Customer No. \_\_\_\_\_ Purchase Order No. \_\_\_\_\_ Fax \_\_\_\_\_

Billing Address \_\_\_\_\_ Patient Name \_\_\_\_\_

Shipping Address \_\_\_\_\_ Order Date \_\_\_\_\_ Measured By \_\_\_\_\_


Telephone \_\_\_\_\_ Email \_\_\_\_\_ Shipping Method \_\_\_\_\_

Credit Card Info \_\_\_\_\_

**arm kit:** Check the box next to the corresponding length and width.

\* If cB is > 30cm hand wrap will not fit, though reduction arm component may still be used.

	Length		Width			Qty.
	Standard	Long	Regular	Wide	X-Wide	
Left	$\ell$ C-G	$\ell$ C-G	cF < 50cm	cF 50-65cm	cF 65-80cm	
Right	< 45cm	> 45cm	cB < 30cm*	cB < 30cm*	cB < 30cm*	



Each arm kit includes:

- 1 reduction arm component
- 1 pair undersleeve arm
- 1 customizable hand wrap
- 1 Built-In-Tension guide card
- 1 direction for use

**Arm Measurement Guide**

G	Axilla
F	Midpoint of E and G
E	Elbow crease
D	Midpoint of C and E
C	Wrist
B	Base of thumb webbing
$\ell$	length in a straight line
c	circumference

**reduction kit glove\*:** Indicate quantity in appropriate check box.

	Small	Medium	Large	X-Large	XX-Large
cB	18-21	20-23	22-25	24-27	26-29
Left					
Right					

\*not included in arm kit

**additional ancillary items:**

undersleeve arm standard Standard <50cm max circ. Quantity (Each)	undersleeve arm wide Wide <80cm max circ. Quantity (Each)	reduction kit shelf strap Quantity (Each)