

## Eclipse®

### Description

Eclipse® is a super absorbent dressing with Crystalock® technology and non-strikethrough backing.

### Indications

Moderate to heavily exuding wounds that require compression or retention bandaging.



How much fluid does Eclipse hold?



Product	SKU	Sizes	Pack Sizes	HCPCS* Code	mL
<b>Eclipse®</b>	CR4503	3in x 3in	10	A6196	69
Super absorbent dressing	CR3818	4in x 4in	20	A6196	122
	CR3769	6in x 6in	20	A6197	275
	CR4439	8in x 8in	10	A6197	488
	CR3743	8in x 12in	20	A6198	732
	CR3808	24in x 16in	10	A6198	2928

## Eclipse®

### ADHERENT

### Description

Eclipse® Adherent is a super absorbent dressing with Crystalock® technology, non-strikethrough backing and soft silicone contact layer.

### Indications

Moderate to heavily exuding wounds that benefit from soft-adherence and require compression or retention bandaging.



How much fluid does Eclipse hold?

Product	SKU	Sizes	Pack Sizes	HCPCS* Code	mL
<b>Eclipse® Adherent</b>	CR4504	3in x 3in	10	A6196	69
Super absorbent dressing with silicone contact layer	CR3881	4in x 4in	10	A6196	122
	CR3883	4in x 8in	10	A6197	244
	CR3863	6in x 6in	10	A6197	275
	CR3864	8in x 12in	10	A6198	732

## Eclipse®

### BORDER

### Description

Eclipse® Border and Eclipse® Border Oval are super absorbent dressings with Crystalock® technology, non-strikethrough backing, soft silicone contact layer and border.

### Indications

Moderate to heavily exuding wounds suitable for a bordered dressing.



How much fluid does Eclipse hold?

Product	SKU	Sizes	Pack Sizes	HCPCS* Code	mL
<b>Eclipse® Border</b>	CR4523	5in x 5in	10	A6196	191
Super absorbent dressing with silicone contact layer and border	CR4296	6in x 6in	10	A6197	275
	CR4297	8in x 12in	10	A6198	732

← New Size Available!

### Eclipse® Border Oval

Super absorbent dressing with silicone contact layer and border

CR4300	4in x 8in	10	A6197	244
CR4299	6in x 8in	10	A6197	366
CR4298	8in x 12in	10	A6197	732

## Eclipse® FOOT

**Description**  
Eclipse® Foot is a super absorbent dressing with Crystalock® technology and non-strikethrough backing. It is uniquely shaped to fit securely around the foot, which is often difficult to dress with regular shaped dressings.

**Indications**  
Moderate to heavily exuding wounds on the foot.



How much fluid does Eclipse hold?

Product	SKU	Sizes	Pack Sizes	HCPCS* Code	mL
<b>Eclipse® Foot</b> Super absorbent foot wrap dressing	CR4219	13in x 19in	5	A6198	1932

## Eclipse® BOOT

**Description**  
Eclipse® Boot is a super absorbent dressing with Crystalock® technology and non-strikethrough backing.

It is uniquely shaped to fit securely around the lower leg.

**Indications**  
Moderate to heavily exuding wounds on the lower leg.



How much fluid does Eclipse hold?

Product	SKU	Sizes	Pack Sizes	HCPCS* Code	mL
<b>Eclipse® Boot</b> Super absorbent leg wrap dressing	CR4237	22in x 19in (SM)	5	A6198	3154
	CR4014	24in x 28in (LG)	5	A6198	5124
	CR4229	28in x 31in (XL)	5	A6198	6930



## Eclipse® CONTOUR

**Description**  
Eclipse® Contour is a super absorbent dressing with Crystalock® technology and non-strikethrough backing. It is uniquely shaped to mold to body contours ensuring optimum contact with the wound.

**Indications**  
Moderate to heavily exuding wounds in difficult to dress areas like the underarm, abdomen, back lower leg and thigh.



How much fluid does Eclipse hold?

Product	SKU	Sizes	Pack Sizes	HCPCS* Code	mL
<b>Eclipse® Contour</b> Super absorbent multisite dressing	CR4398	12in x 20in	10	A6198	1867

## Eclipse® ADHERENT SACRAL

**Description**  
Eclipse® Adherent Sacral is a super absorbent dressing with Crystalock® technology, non-strikethrough backing and soft silicone contact layer.

**Indications**  
Moderate to heavily exuding wounds on the sacral area.



How much fluid does Eclipse hold?

Product	SKU	Sizes	Pack Sizes	HCPCS* Code	mL
<b>Eclipse® Adherent Sacral</b> Super absorbent sacral dressing with silicone contact layer	CR3985	6.5in x 7.5in	10	A6197	394
	CR3986	8.5in x 9in	10	A6198	617

